

> Microscopy Supplies, EasyDip&#8482; Slide Staining Jars

## Microscopy Supplies, EasyDip&#8482; Slide Staining Jars



### Microscopy Supplies, EasyDip&#8482; Slide Staining Jars

The EasyDip&#8482; Slide Staining System has two components: a square staining jar and a 12 - position vertical slide rack. Jars can be loosely joined to each other laterally, therefore making sure they are kept in the same order when moved around on the lab counter. As an extra benefit, they are available in 5 different colors to help better identifying contents or applications. The staining jar is made of resistant acetal plastic and will not break like most glass jars do. It will resist most staining agents including alcohol and xylene (but not phenol, iodine or ferric chloride ). The wide stable base offers greater stability while the inside is recessed, allowing for a smaller reagent volume of only 80 ml. It is easy to clean with no metals to corrode. Ideal for special stains, frozen sections and special processes. Will resist temperatures between -170&deg;C and +121&deg;C. Autoclavable.

Stock #	Description	Price/Cs	Case Pk
M900-12B	EasyDip&#8482; Slide Staining Jars, Blue	\$56.10	6/cs
M900-12G	EasyDip&#8482; Slide Staining Jars, Green	\$56.10	6/cs
M900-12P	EasyDip&#8482; Slide Staining Jars, Pink	\$56.10	6/cs
M900-12W	EasyDip&#8482; Slide Staining Jars, White	\$56.10	6/cs
M900-12Y	EasyDip&#8482; Slide Staining Jars, Yellow	\$56.10	6/cs