

Whirl-Pak Bags, Stand Up Whirl-Pak Sample Bags, Sterile



Whirl-Pak, Stand-Up Whirl-Pak sample bags make sample handling and testing much easier as they give users the ability to perform important laboratory functions with use of both hands -- this plastic sample bag will stand completely by itself, racks or holders are not needed to keep the sample bag upright. When a sample is placed inside, the special gusseted bottom flattens out, allowing the Whirl-Pak bag to stand up on its own.

The white writing area makes sample identification easier and more convenient. The writing strip will accept ordinary ballpoint pen, pencil, or felt-tip marker. The paper-like finish virtually eliminates ink smearing, and provides a surface that will accept a variety of writing instruments.

Capacity is measured in fluid ounces when the stand-up sample bag is closed and tab is folded over three times. Volume and dimensions are approximate; sample bags should not be used at temperatures above 180°F (82°C). Bags can be frozen to any temperature, even liquid nitrogen - careful handling required after freezing. All stand-up Whirl-Pak sample bags are sterilized after manufacturing.

Description	Stock #	Case Pk	Price/Cs
4 oz Whirl-Pak, Stand Up Sample Bags, Sterile	0908-04B	500/cs	\$55.00